



Certification...



certificate can be provided

Both an electronic and hard copy of the BEVA

TKP qPCR (triple test) Taylorella, Klebsiella & Pseudomonas

TKP q PCR a rapid assay with **same day** results using QIAGEN state of the art technology giving accurate sensitive results for all 3 CEM organisms. For PCR swabs can be plain or in charcoal medium, Same day results for samples received by 10.00 am



Aerobic and microaerophilic

Tradition 7 day microaerophilic and aerobic culture to examine for all 3 CEM organisms Amies charcoal swabs only which must arrive at the laboratory within 48 hours of collection, please ensure swabs are within their expiry date.

All CEM tests are accredited and quality assessed by the VETQAS proficiency scheme.



Mare endometrial swabs gel or plain swabs

Cytology report- microscopic smear, including examination for the presence of PMNs, and bacteria, same day report

Overnight aerobic culture for Klebsiella & Pseudomonas, report at 24 & 48hrs

Serology

EVA & EIA both can be performed by ELISA, Enzyme Linked Immunosorbent Assay (not for export), results are available in 48 hours

For export: EIA a Coggins blood test is required,
EVA is sent to the Animal & Plant Health Agency, Weybridge

Vaccinated Stallion and teasers require a VNT test, results available in 5 days
4ml clotted blood / 2ml serum

Further information and swabbing techniques can be found at:

www.codes.hblb.org.uk and www.beva.org.uk

For supplies of Pre-paid postage labels and custom printed submission forms please email your requirements to info@jsc-equine-laboratory.co.uk